

Product datasheet for **AP54800PU-N**

BPGM Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/50 - 1/500. ELISA. Western Blot: 1/200 - 1/2000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of human BPGM protein
Specificity:	This antibody reacts with 30 kDa BPGM protein.
Formulation:	Tris 0,1M, glycine 0,1M, sucrose 2% State: Purified State: Lyophilized purified powder
Reconstitution Method:	Restore in distilled water.
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	bisphosphoglycerate mutase
Database Link:	Entrez Gene 669 Human P07738
Synonyms:	BPGM



[View online »](#)