

Product datasheet for **AP54802SU-N**

BRCA2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/50 - 1/500. ELISA. Western Blot: 1/200 - 1/2000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from N-term domain of human BRCA2
Specificity:	This antibody reacts with 384 kDa BRCA2 protein. Cross reacts with phospho-Ser 28 BRCA2.
Formulation:	State: Serum State: Lyophilized serum Preservative: None
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	BRCA2, DNA repair associated
Database Link:	Entrez Gene 675 Human P51587
Synonyms:	FACD



[View online »](#)