

Product datasheet for **AP54804SU-N**

Btn2a2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, WB
Recommended Dilution:	ELISA. Western Blot. Immunohistochemistry.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptides derived from internal domain of Mouse Btn2a2 protein
Specificity:	This antibody reacts with 55 kDa Btn2a2 protein.
Formulation:	State: Serum State: Lyophilized serum Preservative: None
Reconstitution Method:	Restore in distilled water. Recommended: 0.2 ml.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	butyrophilin, subfamily 2, member A2
Database Link:	Entrez Gene 238555 Mouse A4QPC6
Synonyms:	Butyrophilin subfamily 2 member A2, BT2.2, BTF2



[View online »](#)