

## Product datasheet for **AP54814SU-N**

### adhA Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	<b>ELISA.</b> <b>Western Blot:</b> 1/500 - 1/2000.
Reactivity:	Burkholderia
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from N-ter domain of the Burkholderia adhA protein
Specificity:	This antibody reacts with 22 kDa adhA protein of different strains of Burkholderia.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	<a href="#">Q3JS59</a>



[View online »](#)