

Product datasheet for **AP54818SU-N**

CLSTN1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot: 1/500-1/5000. Immunohistochemistry: 1/200-1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of CSTN1
Specificity:	This antibody reacts with 109 kDa CSTN1 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	calsyntenin 1
Database Link:	Entrez Gene 22883 Human O94985
Synonyms:	CLSTN1, CS1, KIAA0911, Alcadin-alpha, Cadherin XB31 alpha, Alc-alpha



[View online »](#)