

Product datasheet for **AP54819SU-N**

TPTE Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	ELISA. Western Blot. Immunofluorescence. Immunohistochemistry on Paraffin Sections.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptides derived from internal domains of the Human TPTE protein
Specificity:	Reacts with Human 64 kDa TPTE protein. May cross react with mouse and other species due to sequence homology.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	transmembrane phosphatase with tensin homology
Database Link:	Entrez Gene 7179 Human P56180
Synonyms:	BJ-HCC-5



[View online »](#)