

## Product datasheet for **AP54821SU-N**

### CPA4 Rabbit Polyclonal Antibody

#### Product data:

|                        |   |
|------------------------|---|
| Product Type:          | Primary Antibodies  |
| Applications:          | ELISA, IHC, WB  |
| Recommended Dilution:  | <b>Immunohistochemistry:</b> 1/100 - 1/500.<br><b>ELISA.</b><br><b>Western Blot:</b> 1/200 - 1/1000.  |
| Reactivity:            | Human   |
| Host:                  | Rabbit  |
| Clonality:             | Polyclonal  |
| Immunogen:             | Synthetic peptide derived from an internal domain of Human CPA4   |
| Specificity:           | This antibody reacts with 47KDa CPA4.   |
| Formulation:           | State: Serum<br>State: Lyophilized serum  |
| Reconstitution Method: | Restore in distilled water.   |
| Conjugation:           | Unconjugated  |
| Storage:               | Prior to reconstitution store at -20°C.<br>Following reconstitution store undiluted at 2-8°C for one month<br>or (in aliquots) at -20°C for longer.<br>Avoid repeated freezing and thawing. |
| Stability:             | Shelf life: one year from despatch.   |
| Gene Name:             | carboxypeptidase A4   |
| Database Link:         | <a href="#">Entrez Gene 51200 Human Q9UI42</a>  |
| Synonyms:              | CPA4, CPA3  |



[View online »](#)