

Product datasheet for **AP54823SU-N**

CMYA5 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/50 – 1/500. ELISA. Western Blot: 1/200 to 1/2000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of human CMYA5
Specificity:	This antibody reacts with 36 kDa CMYA5 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	cardiomyopathy associated 5
Database Link:	Entrez Gene 202333 Human Q8N3K9
Synonyms:	Dystrobrevin-binding protein 2, Genethonin-3, CMYA5, C5orf10, DTNBP2, SPRYD2, TRIM76

