

## Product datasheet for **AP54854SU-N**

### Collagen IV (COL4A4) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>Immunohistochemistry:</b> 1/50 - 1/500. <b>ELISA.</b> <b>Western Blot:</b> 1/500 - 1/2000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of human COL4A4
Specificity:	This antibody reacts with precursor 164kDa COL4A4 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	collagen type IV alpha 4 chain
Database Link:	<a href="#">Entrez Gene 1286 Human P53420</a>
Synonyms:	COL4A4, Collagen alpha-4(IV) chain



[View online »](#)