

Product datasheet for **AP54881SU-N**

Prkdc Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/50 - 1/500. ELISA. Western Blot: 1/100 - 1/1000.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from internal domain of mouse PRKDC protein
Specificity:	This antibody reacts with 471 kDaPRKDC protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	protein kinase, DNA activated, catalytic polypeptide
Database Link:	Entrez Gene 19090 Mouse P97313
Synonyms:	DNPK1, p460, HYRC, HYRC1



[View online »](#)