

Product datasheet for **AP54903SU-N**

Dytn Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot: 1/500-1/5000. Immunocytochemistry: 1/100-1/500.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from internal domain of the Mouse Dystrotelin protein
Specificity:	Reacts with Mouse 74 kDa Dystrotelin (Dytn) protein.
Formulation:	State: Serum State: Lyophilized powder Preservative: None
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	dystrotelin
Database Link:	Entrez Gene 241073 Mouse A2C198
Synonyms:	DYTN



[View online »](#)