

Product datasheet for **AP54910SU-N**

ERC1 Rabbit Polyclonal Antibody

Product data:

| | |
|------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | ELISA, IHC, WB |
| Recommended Dilution: | ELISA. Immunohistochemistry: 1/50 - 1/500. Western Blot: 1/200 - 1/2000. |
| Reactivity: | Human, Mouse |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | Synthetic peptide derived from Nter domain of ELKS protein |
| Specificity: | This antibody reacts with human and mouse ERC-1. |
| Formulation: | State: Serum State: Lyophilized Serum |
| Reconstitution Method: | Restore in distilled water. |
| Conjugation: | Unconjugated |
| Storage: | Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: one year from despatch. |
| Gene Name: | ELKS/RAB6-interacting/CAST family member 1 |
| Database Link: | Entrez Gene 23085 Human Q8IUD2 |
| Synonyms: | ERC protein 1, KIAA1081, RAB6IP2 |



[View online »](#)