

Product datasheet for **AP54912SU-N**

ETFDH Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Immunohistochemistry: 1/50 - 1/500. Western Blot: 1/200 - 1/2000.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Nter domain of human ETFD
Specificity:	This antibody reacts with 68 kDa ETFDH protein.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	electron transfer flavoprotein dehydrogenase
Database Link:	Entrez Gene 2110 Human Q16134
Synonyms:	ETF-ubiquinone oxidoreductase, ETF-QO



[View online »](#)