

Product datasheet for **AP54919SU-N**

Ephrin B3 (EFNB3) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Immunohistochemistry on frozen sections: 1/50 - 1/500. Western Blot: 1/500 - 1/2000.
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Nter domain of human EFNB3
Specificity:	This antibody reacts with 37 kDa EFNB3 protein.
Formulation:	State: Serum State: Lyophilized Serum
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	ephrin B3
Database Link:	Entrez Gene 1949 Human Q15768
Synonyms:	LERK-8, ELK-L3, EPLG8, LERK8, EFNB3



[View online »](#)