

## Product datasheet for **AP54920SU-N**

### EMP3 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunohistochemistry on frozen sections:</b> 1/100 - 1/1000. <b>Western Blot:</b> 1/1000 - 1/5000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from human EMP3 protein
Specificity:	This antibody reacts with 18 kDa EMP3 protein.
Formulation:	State: Serum State: Lyophilized Serum
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	epithelial membrane protein 3
Database Link:	<a href="#">Entrez Gene 2014 Human P54852</a>
Synonyms:	Epithelial membrane protein 3, HNMP-1



[View online »](#)