


Product datasheet for **AP54921SU-N**

EPB42 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot: 1/200-1/1000. Immunohistochemistry on Frozen Sections : 1/200-1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A 15 residue synthetic peptide derived from specific domain of band EP4.2 protein
Specificity:	This antibody reacts with P4.2 from erythrocytes. Does not cross react with transglutaminases
Formulation:	State: Serum State: Lyophilized Serum Preservative: None
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	erythrocyte membrane protein band 4.2
Database Link:	Entrez Gene 2038 Human P16452
Synonyms:	Band 4.2, Erythrocyte protein 4.2, P4.2



Product images:

Western Blot of Human Erythrocyte extract
using EPB42 / E42P Antibody Cat.-No
AP54921SU-N