

Product datasheet for **AP54938SU-N**

Flagellin Fla1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot: 1/500-1/2000. Immunohistochemistry.
Reactivity:	Roseburia hominis
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from internal domain of Flagellin-1 protein
Specificity:	This antibody reacts with 29 kDa Fla1 from Roseburia hominis (clostridia).
Formulation:	State: Serum State: Lyophilized Serum Preservative: None
Reconstitution Method:	Restore in 200 µl distilled water.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: 6 month from despatch.
Database Link:	A7IZH9



[View online »](#)