

Product datasheet for **AP54939SU-N**

Flagellin Fla2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Immunohistochemistry. Western Blot: 1/500-1/2000.
Reactivity:	Roseburia hominis
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from an internal domain of Flagellin 2 protein
Specificity:	This antibody reacts with 31 kDa Flagellin 2 from Roseburia hominis (clostridia).
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in 0.2 ml distilled water.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	A7IZIO



[View online »](#)