

## Product datasheet for **AP54940PU-N**

### FSTL1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunohistochemistry:</b> 1/50 - 1/500. <b>Western Blot:</b> 1/500 to 1/2000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Nter domain of FSTL1 protein
Specificity:	This antibody reacts with 35 kDa FSTL1 protein. Cross reacts with mouse, rat and other species due to sequence homology.
Formulation:	Tris 0,1M, glycine 0,1M, sucrose 2% State: Purified State: Lyophilized purified powder
Reconstitution Method:	Restore in distilled water.
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	follistatin like 1
Database Link:	<a href="#">Entrez Gene 11167 Human Q12841</a>
Synonyms:	FSTL1, FRP, Follistatin-like protein 1



[View online »](#)