

## Product datasheet for **AP54941SU-N**

### FOXD3 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	<b>ELISA.</b> <b>Western Blot:</b> 1/500 - 1/5000.
Reactivity:	Chicken
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from N-term domain of the FOXD3 chicken protein
Specificity:	This antibody reacts with 40 kDa FOXD3 protein. Cross reacts with protein from mouse, human and other species due to sequence homology.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	forkhead box D3
Database Link:	<a href="#">Entrez Gene 395794 Chicken P79772</a>
Synonyms:	Forkhead box protein D3



[View online »](#)