

## Product datasheet for **AP54946SU-N**

### **FXR1 Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunohistochemistry:</b> 1/100 - 1/1000. <b>Western Blot:</b> 1/1000 - 1/5000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Nter domain of human FXR1 protein
Specificity:	This antibody reacts with 70 kDa hFXR1p protein. Cross reacts with protein from mouse, rat and other species due to sequence homology.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	FMR1 autosomal homolog 1
Database Link:	<a href="#">Entrez Gene 8087 Human P51114</a>
Synonyms:	Fragile X mental retardation syndrome-related protein 1, hFXR1p



[View online »](#)