

Product datasheet for **AP54947SU-N**

Fli1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot: 1/1000 - 1/5000. Immunohistochemistry: 1/100 - 1/1000.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from internal domain of mouse FLI-1 protein
Specificity:	This antibody reacts with 51 kDa FLI-1 protein.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	Friend leukemia integration 1
Database Link:	Entrez Gene 14247 Mouse P26323
Synonyms:	Fli-1 proto-oncogene, ERGB transcription factor, vascular tumor marker



[View online »](#)