

Product datasheet for **AP54953PU-N**

Slc6a1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Immunohistochemistry: 1/5 - 1/50. Western Blot: 1/20 - 1/500.
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A 20 residue synthetic peptide derived from C-terminal domain of rat GABA transporter-1
Specificity:	This antibody reacts with 72 kDa GABA transporter-1.
Formulation:	Tris 0.1 M, glycine 0.1 M, Glycerol 50%, (pH 7.8) State: Aff - Purified State: Liquid purified Ig fraction Preservative: 0.02% as preservative
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	solute carrier family 6 member 1
Database Link:	Entrez Gene 79212 Rat P23978
Synonyms:	GABATR, GABT1, GAT-1, GABA transporter 1



[View online »](#)