

Product datasheet for **AP54959SU-N**

Rtl3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blot.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptides derived from internal domain of the mouse Zcchc5 protein.
Specificity:	This antibody reacts with 61 kDa Zcchc5 protein.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	zinc finger, CCHC domain containing 5
Database Link:	Entrez Gene 213436 Mouse Q6P1Y1
Synonyms:	FLJ38865; Mar3; Mart3; ZHC5



[View online »](#)