

Product datasheet for **AP54967PU-N**

Gastrin (GAST) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, WB
Recommended Dilution:	ELISA. Western Blot. Immunohistochemistry.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptides derived from internal domain of the human Gastrin precursor protein
Specificity:	This antibody reacts with 11 kDa GAST protein.
Formulation:	Tris 0,1M, glycine 0,1M, sucrose 2% State: Purified State: Lyophilized purified powder
Reconstitution Method:	Restore in distilled water.
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	gastrin
Database Link:	Entrez Gene 2520 Human P01350
Synonyms:	GAST, GAS



[View online »](#)