

Product datasheet for **AP54971SU-N**

GLUT4 (SLC2A4) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/200 - 1/1000. ELISA. Western Blot: 1/500 - 1/5000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of human GluT4
Specificity:	This antibody reacts with 54 kDa GTR4 protein.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	solute carrier family 2 member 4
Database Link:	Entrez Gene 6517 Human P14672
Synonyms:	GLUT-4, Glucose transporter 4, Glucose transporter type 4 insulin-responsive



[View online »](#)