

Product datasheet for **AP54974SU-N**

Gpr119 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	ELISA. Western Blot: 1/1000 - 1/5000. Immunohistochemistry: 1/100 - 1/1000.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain mouse Gpr119 protein
Specificity:	This antibody reacts with Gpr 119 protein.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	G-protein coupled receptor 119
Database Link:	Entrez Gene 236781 Mouse Q7TQP3
Synonyms:	G-protein coupled receptor 119



[View online »](#)