

## Product datasheet for **AP54980SU-N**

### Glutathione Peroxidase 2 (GPX2) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunohistochemistry:</b> 1/100 - 1/1000. <b>Western Blot:</b> 1/250 - 1/5000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of GPx-2 protein
Specificity:	This antibody reacts with 22 kDa GPx-2 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	glutathione peroxidase 2
Database Link:	<a href="#">Entrez Gene 2877 Human P18283</a>
Synonyms:	GPx-2, GSHPx-2, GPRP, Gastrointestinal glutathione peroxidase



[View online »](#)