

Product datasheet for **AP54997SU-N**

Hist1h1a Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blot: 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from N-term domain of human histone H1.1
Specificity:	This antibody reacts with recombinant H1.1. Does not cross react with other human H1 isoforms or core histones.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	histone cluster 1, H1a
Database Link:	P43275
Synonyms:	HIST1H1A, H1F1



[View online »](#)

Product images:



Silver staining
recombinant
human H1

Western Blot of Histone H1.1. antibody using AP54997SU-N at 1/2000 dilution



H1.1 Isoform