

Product datasheet for **AP54999SU-N**

Histone H1.3 (HIST1H1D) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, WB
Recommended Dilution:	Immunocytochemistry: 1/50. Immunofluorescence: 1/50. ELISA. Western Blot: 1/500.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from N-term domain of human histone H1.3 protein.
Specificity:	This antibody reacts with H1.3.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	histone cluster 1, H1d
Database Link:	Entrez Gene 3007 Human P16402
Synonyms:	HIST1H1D, H1F3, Histone H1c



[View online »](#)