

Product datasheet for **AP55001SU-N**

H1oo (H1FOO) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot: 1/200-1/2000. Immunohistochemistry: 1/50-1/500.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from the C-terminal domain of H1oo
Specificity:	This antibody reacts with Mouse 32 kDa Histone H1oo.
Formulation:	State: Serum State: Lyophilized Serum Preservative: None
Reconstitution Method:	Restore in 0.2 ml distilled water.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	H1 histone family member O, oocyte specific
Database Link:	Q8IZA3
Synonyms:	H1OO, OSH1, Oocyte-specific histone H1



[View online »](#)