

Product datasheet for **AP55002SU-N**

H1oo (H1FOO) pSer20 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blot: 1/200 - 1/2000.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Syntheticpeptide derived from N-term domain of H100 containing pS20
Specificity:	This antibody reacts with 32 kDa Histone H1oo phosphorylated at serine 20.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	H1 histone family member O, oocyte specific
Database Link:	Q8IZA3
Synonyms:	H100, OSH1, Oocyte-specific histone H1



[View online »](#)