

Product datasheet for **AP55003PU-N**

Histone H2A.X (H2AFX) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	ELISA. Western Blot. Immunofluorescence. Immunohistochemistry.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C term part of the human H2a/x protein
Specificity:	This antibody reacts with the 15kDa H2a/x protein. May cross react with H2A phosphoserine-140 and with bovine, mouse and rat proteins.
Formulation:	Tris 0,1M, glycine 0,1M, sucrose 2% State: Purified State: Lyophilized purified powder
Reconstitution Method:	Restore in distilled water.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	H2A histone family member X
Database Link:	Entrez Gene 3014 Human P16104
Synonyms:	H2a/x, H2AFX, H2AX



[View online »](#)