

Product datasheet for **AP55018SU-N**

Hyaluronan synthase Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blot: 1/200 - 1/2000.
Reactivity:	Paramecium Bursaria Chlorella Virus 1
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of Hyaluronan synthase from Paramecium bursaria chlorella virus-1 (PBCV-1)
Specificity:	This antibody reacts with 65 kDa Hyaluronan synthase from PBCV-1.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	Q9WFS9



[View online »](#)