

## Product datasheet for **AP55020SU-N**

### **PAMR1 Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Western Blot.</b> <b>Immunohistochemistry.</b>
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from internal domain of human PAMR1 protein
Specificity:	This antibody reacts with 80 kDa PAMR1 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	peptidase domain containing associated with muscle regeneration 1
Database Link:	<a href="#">Entrez Gene 25891 Human</a> <a href="#">Q6UXH9</a>
Synonyms:	RAMP, DKFZP586H2123, FP938, UNQ699/PRO1344



[View online »](#)