

## Product datasheet for **AP55022SU-N**

### Pyrophosphatase 1 (PPA1) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>Immunohistochemistry:</b> 1/50 - 1/200. <b>ELISA.</b> <b>Western Blot:</b> 1/200 - 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of PPA1 protein
Specificity:	This antibody reacts with 32 kDa PPA1 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	pyrophosphatase (inorganic) 1
Database Link:	<a href="#">Entrez Gene 5464 Human</a> <a href="#">Q15181</a>
Synonyms:	Inorganic pyrophosphatase, Pyrophosphate phospho-hydrolase, PPase, IOPPP, PP



[View online »](#)