

Product datasheet for **AP55044SU-N**

Katnal1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/50 - 1/200. ELISA. Western Blot: 1/100 - 1/1000.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of human KATL1 protein
Specificity:	This antibody reacts with 55 kDa KATL1 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	katanin p60 subunit A-like 1
Database Link:	Entrez Gene 231912 Mouse Q8K0T4
Synonyms:	p60 katanin-like 1



[View online »](#)