

Product datasheet for **AP55054SU-N**

LLGL2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/100 - 1/1000. ELISA. Western Blot: 1/1000 - 1/5000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Nter domain of human L2GL2 protein
Specificity:	This antibody reacts with 113 kDa L2GL2 protein.
Formulation:	0.1M Tris, 0.1M Glycine, 2% Sucrose State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	LLGL2, scribble cell polarity complex component
Database Link:	Entrez Gene 3993 Human Q6P1M3
Synonyms:	HGL, Lethal(2) giant larvae protein homolog 2



[View online »](#)