

Product datasheet for **AP55059SU-N**

LGR6 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot: 1/200 - 1/1000. Immunohistochemistry: 1/50 - 1/500.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from human LGR6 protein
Specificity:	This antibody reacts with the LGR6 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	leucine-rich repeat containing G protein-coupled receptor 6
Database Link:	Entrez Gene 59352 Human Q9HBX8
Synonyms:	G-protein coupled receptor 6



[View online »](#)