

Product datasheet for **AP55064SU-N**

Lysophospholipase 1 (LYPLA1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/200 - 1/1000. ELISA. Western Blot: 1/500 - 1/5000.
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of Lyso PLA 1 protein
Specificity:	This antibody reacts with 24 kDa LYPA1 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	lysophospholipase I
Database Link:	O75608
Synonyms:	Acyl-protein thioesterase 1, APT1, APT-1, LPL1, LPL-1, Lysophospholipase 1, LysoPLA I



[View online »](#)