

Product datasheet for **AP55065SU-N**

Marveld3 (Isoform 2) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Immunohistochemistry. Western Blot: 1/200 - 1/1000.
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of mouse MALD3 isoform-2 protein
Specificity:	This antibody reacts with 43 kDa MALD3 isoform-2 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	MARVEL (membrane-associating) domain containing 3
Database Link:	Entrez Gene 73608 Mouse Q9D956
Synonyms:	MARVEL domain-containing protein 3, MRVLDC3



[View online »](#)