

Product datasheet for **AP55066SU-N**

Marveld3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/50 - 1/500. ELISA. Western Blot: 1/500 - 1/2000.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of MALD3 protein
Specificity:	This antibody reacts with 42 kDa MALD3 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	MARVEL (membrane-associating) domain containing 3
Database Link:	Entrez Gene 73608 Mouse Q9D956
Synonyms:	MARVEL domain-containing protein 3, MRVLDC3



[View online »](#)