

## Product datasheet for **AP55069SU-N**

### **MKRN1 Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>Immunohistochemistry:</b> 1/100 - 1/1000. <b>ELISA.</b> <b>Western Blot:</b> 1/500 - 1/2000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from N-term domain of human MKRN1 protein
Specificity:	This antibody reacts with 53 kDa MKRN1 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	makorin ring finger protein 1
Database Link:	<a href="#">Entrez Gene 23608 Human Q9UHC7</a>
Synonyms:	MKRN1, RNF61, RING finger protein 61



[View online »](#)