

Product datasheet for **AP55099SU-N**

NOX3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/50 - 1/500. ELISA. Western Blot: 1/200 - 1/2000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from an internal domain of human NOX3
Specificity:	This antibody reacts with 65 kDa NOX3 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	NADPH oxidase 3
Database Link:	<u>Entrez Gene 50508 Human Q9HBY0</u>
Synonyms:	MOX2, Mitogenic oxidase 2, MOX-2, gp91phox homolog 3, GP91-3



[View online »](#)