

Product datasheet for **AP55101SU-N**

UNC5A Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, WB
Recommended Dilution:	ELISA. Western Blot. Immunohistochemistry.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptides derived from internal domain of the human UNC5A protein.
Specificity:	This antibody reacts with 93kDa UNC5A protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	unc-5 netrin receptor A
Database Link:	Entrez Gene 90249 Human Q6ZN44
Synonyms:	Protein unc-5 homolog A, Unc-5 homolog 1, KIAA1976, UNC5H1



[View online »](#)