

## Product datasheet for **AP55102SU-N**

### PPP1R9A Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>Immunohistochemistry:</b> 1/50 - 1/500. <b>ELISA.</b> <b>Western Blot:</b> 1/200 - 1/2000.
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from N-term domain of rat NEB1
Specificity:	This antibody reacts with 122 kDa NEB1 protein.
Formulation:	State: Serum State: Lyophilized serum Preservative: None
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	protein phosphatase 1, regulatory subunit 9A
Database Link:	<a href="#">Entrez Gene 84685 Rat O35867</a>
Synonyms:	FLJ20068; KIAA1222; Neurabin-I; NRB1; NRBI; OTTHUMP00000208240



[View online »](#)