

## Product datasheet for **AP55116SU-N**

### NUCKS1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, WB
Recommended Dilution:	<b>ELISA.</b> <b>Western Blot.</b> <b>Immunohistochemistry.</b>
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptides derived from internal domain of the NUCKS1
Specificity:	This antibody reacts with 27 kDa NUCKS1 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	nuclear casein kinase and cyclin-dependent kinase substrate 1
Database Link:	<a href="#">Entrez Gene 64710 Human Q9H1E3</a>
Synonyms:	NUCKS1, JC7, P1



[View online »](#)