

Product datasheet for **AP55119SU-N**

Okadaic Acid Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1/5000 - 1/20000.
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Ovalbumin-Okadaic Acid
Specificity:	This antibody reacts with Okadaic acid from Halichondria okadai.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)