

Product datasheet for **AP55124SU-N**

Orexin Receptor 2 (HCRT2) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/50 - 1/200. Western Blot: 1/100 - 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of human Ox2 receptor protein
Specificity:	This antibody reacts with 50 kDa OX2 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	hypocretin receptor 2
Database Link:	Entrez Gene 3062 Human O43614
Synonyms:	HCRT2, Ox2r, Hypocretin receptor type 2



[View online »](#)