

## Product datasheet for **AP55128SU-N**

### PAR6 (C-term) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	<b>ELISA.</b> <b>Western Blot:</b> 1/200 - 1/1000.
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C terminal region of PAR6 protein
Specificity:	This antibody reacts with 45 kDa PAR6 (Cter Domain) of Phallusia Mammilata.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	<a href="#">Q52PI6</a>



[View online »](#)